alber®

erivo®

Gebrauchsanweisung Erivo carbon

DE

Instructions for Use Erivo carbon

EN

CE | UK
CA

1596329

Komponenten

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Components


- | | | | |
|---|--------------|----|----------------------|
| 1 | Headlight | 8 | Parking brake |
| 2 | Control unit | 9 | Anti-tip wheels |
| 3 | Armrest | 10 | Drive wheel |
| 4 | Backrest | 11 | Storage basket |
| 5 | Push handle | 12 | Caster (front wheel) |
| 6 | Bag | 13 | Footrest |
| 7 | Battery pack | 14 | Seat cushion |

Charger

- 50 Charger
- 51 LED indicator
- 52 Charging plug
- 53 Mains plug





Battery Pack

- 7 Battery pack
- 49 Charging socket
- 61 System power socket

Control Unit

- 20 On/Off button / Anti-theft lock
- 21 Emergency call button
- 22 Horn button
- 23 Screen (*only visible when the device is switched on*)
- 24 Joystick
- 25 Drive programme selector button
- 26 Button for selecting a lower speed level
- 27 Button for selecting a higher speed level

47.0001.4.02.01
 Status: 2026-03-19

This instructions for use (IFU) provides information about the product *Erivo carbon*, hereinafter referred to as *Erivo*.



The current version of these instructions for use is available for download from our website www.alber.de. If you require a version in larger print, please download the IFU as a PDF file from our website and increase the display size on your screen.

The illustrations contained in these instructions for use may differ from the product actually delivered to you.

Read these instructions for use carefully before commissioning the Erivo. They contain all the necessary information and safety instructions for operating the device.

Service Center (Deutschland)

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8.00 - 18.00 Uhr

Freitags von

8.00 - 16.00 Uhr

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Telefon (0800) 9096-250

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Table of Contents

1. Introduction	2	5.6 Triggering an Emergency Call	21
1.1 Intended Purpose	2	5.7 Limiting Maximum Speed (Speed Levels)	22
1.2 Contraindications	2	5.8 Driving in Beginner Mode / Exiting Beginner Mode	23
1.3 Important Safety Instructions – Please observe without exception	2	5.9 Driving with Speed Control	24
1.4 Intended Use of the Erivo	2	5.10 Pushing the Erivo (Push Mode)	25
1.5 Handling Instructions	2	5.11 Releasing / Activating the Parking Brakes	26
1.6 Symbol Explanation	3	5.12 Activating / Deactivating the Anti-Theft Lock	26
1.7 Permitted Operating Conditions / Locations	4	5.13 Screen Displays	27
1.8 Standard Scope of Delivery	5	5.14 Using the Lighting	28
1.9 Optional Accessories	5	5.15 Using the USB Charging Ports	29
2. Commissioning	5	6. Faults	29
2.1 Adapting the Erivo to the User	6	6.1 Fault Displays on the Screen	29
2.1.1 Swapping the Armrests	6	6.2 Other Faults	30
2.1.2 Extending / Shortening the Armrests	6	7. Automatic Self-Shutdown	31
2.2 Fitting / Adjusting the Safety Belt	7	7.1 Self-Shutdown Due to Internal Faults	31
2.3 Attaching / Removing the Bag	7	7.2 Self-Shutdown During Non-Use / Cancelling Self-Shutdown	31
2.4 Using the Storage Basket	8	8. Checking Tyre Pressure	31
2.5 Switching the System On	8	9. Transferring	32
2.6 Switching the System Off	9	10. Removing / Placing the Seat Cushion	32
2.7 Range	9	11. Transport	32
3. Safety and Hazard Notices for Driving	9	11.1 Folding and Unfolding the Wheelchair	32
3.1 General Notes	9	11.1.1 Folding the Wheelchair	33
3.2 Driving Training Notes	9	11.1.2 Unfolding the Wheelchair	33
3.3 Safety Instructions	10	11.2 Safety Instructions for Loading the Wheelchair	34
3.4 Special Safety Instructions for Slopes, Descents and Ramps	11	11.2.1 Loading Without a Ramp or Lifting Device - Folded	34
3.5 Obstacles	12	11.2.2 Loading via Ramp or Lifting Device - Unfolded - by Accompanying Person	34
3.6 Hazardous Locations and Situations	12	11.2.3 Loading via Ramp or Lifting Device - Unfolded - by Wheelchair User	34
3.7 Electromagnetic Interference (EMI)	13	11.3 Transporting the Wheelchair in a Vehicle	34
4. Battery Pack	13	12. Storage, Care, Recommissioning, Maintenance and Disposal	36
4.1 Removing the Battery Pack	14	12.1 Storing the Wheelchair	36
4.2 Fitting the Battery Pack	14	12.2 Caring for the Wheelchair	36
4.3 Charging the Battery Pack	15	12.3 Recommissioning	37
4.3.1 Charging at the Control Unit	15	12.4 Maintenance	37
4.3.2 Charging at the Battery Pack Charging Socket	14	12.5 Disposal	37
4.4 Capacity Indicator	16	13. Product Safety Notices	38
4.5 General Charging Notes	16	14. Product Lifespan	38
4.6 Safety Instructions for the Charger and Charging Process	17	15. Warranty, Guarantee and Liability	38
4.7 Safety Instructions for the Battery Pack	18	15.1 Defect Warranty	38
4.8 Storing the Battery Pack	18	15.2 Durability Guarantee	38
4.9 Safety and Warning Notices for Transport and Shipping of the Battery Pack	18	15.3 Liability	38
4.10 Important Information for Air Travel	19	16. Important Legal Notice for Users of this Product	38
5. Driving the Erivo	19	17. Technical Data	39
5.1 Switching the System On	19	18. Labels and Markings	40
5.2 Switching the System Off	19	Appendix	
5.3 Driving with the Joystick	20	A - INVACARE / Alber Representatives	43
5.4 Using the Horn	21	B - Inspection Checks	43
5.5 Emergency Braking	21		

1. Introduction

1.1 Intended Purpose

The Erivo is an electric wheelchair. It is intended as a medical aid for disabled persons in order to increase their mobility and flexibility. It is further intended for disabilities that do not require individual adjustment of the seating geometry.

1.2 Contraindications

For safety reasons, the Erivo must not be operated by persons whose psychological conditions — whether temporary or permanent — impair their attention or judgement.

1.3 Important Safety Instructions – Please observe without exception

For safety reasons, the Erivo may only be operated by persons who:

- have been instructed in its handling and/or have at least read and understood these complete instructions for use,
- are able to move and coordinate both hands and arms without significant restriction,
- are physically and mentally capable of operating the wheelchair safely in all operating situations (e.g. road traffic) and, in the event of failure of the electrically driven wheels, of braking the wheelchair and bringing it to a safe standstill,
- are not under the influence of medications (which impair reaction times), alcohol, or drugs.

If you are not yet confident in handling the Erivo, please also contact your specialist dealer.

During operation of the Erivo, the prescribed values (e.g. maximum gradient, tyre pressure, etc.) and general operating instructions must be strictly observed. Limit values must under no circumstances be exceeded.

Travel near strong electrical interference fields should be avoided.

In rare cases, operation of the Erivo may possibly interfere with other installations, for example anti-theft systems in department stores. Escalators and travelators must not be used with the Erivo.

During sporting activities performed in the wheelchair, such as weightlifting or similar, the Erivo must be switched off.

It is also not permitted to combine the Erivo with accessories that have not been approved by Alber.

Commissioning or operating the Erivo without reading and understanding these instructions for use is not permitted and constitutes improper use.

In the event of improper use, there is a risk of personal injury.

1.4 Intended Use of the Erivo

The Erivo is intended exclusively for the transportation of persons with mobility impairments. Use and transport, as well as maintenance and service of the Erivo, must be carried out in accordance with the information in these instructions for use.

1.5 Handling Instructions

No parts other than accessories approved by Alber for use may be fitted. Likewise, the Erivo and its accessories must not be technically modified.

The Erivo must be handled under the following conditions:

- Compliance with the information, instructions and recommendations in these instructions for use
- The Erivo is operated exclusively by an instructed person or a person who has read and understood the complete instructions for use and has thus qualified for use.
- No technical modifications have been made to the Erivo by the user or by third parties.

An instructed person is defined as someone who has demonstrably been informed of the tasks assigned to them, the possible hazards in the event of improper behaviour, and has been introduced to the handling of the Erivo, or who has qualified for use by reading and understanding the complete instructions for use.

Operation of the Erivo by uninstructed or unqualified users is expressly prohibited.

The Erivo must not be used for purposes contrary to its intended use. This applies in particular to all types of goods transport, such as the transport of objects or additional persons (e.g. children). Intended use also includes compliance with the maintenance specifications in these instructions for use, and observance of the safety and hazard notices for driving operations.

Alber GmbH regards the following as misuse of the Erivo:

- Use of the device contrary to the instructions and recommendations of these instructions for use
- Exceeding the technical limits defined in these instructions for use
- Technical modifications to or within the device
- Fitting and use of third-party parts not manufactured or offered for use by Alber

For damage resulting from:

- misuse of the device,
- operation by a person not instructed or qualified in handling the device,
- use contrary to the instructions and recommendations of these instructions for use,
- exceeding the technical limits defined in these instructions for use,

Alber GmbH disclaims all liability.

Alber GmbH excludes liability for damages caused by:

- **improper use of the Erivo and its accessories,**
- **handling of the Erivo and its accessories by an uninstructed or unqualified person,**
- **use of the Erivo and its accessories contrary to the instructions, in particular the safety and hazard notices in these instructions for use,**
- **exceeding the technical limits defined in these instructions for use, causing harm to the user or third parties.**

Before commissioning the Erivo, familiarise yourself with the safety and hazard notices in the individual chapters of these instructions for use.

1.6 Symbol Explanation

Important tips and notices within these instructions for use are marked as follows:

Note on tips and special information.

Warning of possible hazards to your safety and health, as well as notices of possible injury risks. Warning of possible technical problems or damage.

Observe these notices and warnings without fail in order to prevent injury to persons and damage to the product!

Information in these instructions for use, for example „front“, „rear“, „left“, „right“, „forward“, „backward“ etc., refers to the position as seen **from the driver's perspective**.

The following are explanations of the symbols used on the labels and in part within these instructions for use.

The Erivo fulfils all applicable chapters of EN 12184 for electrically powered wheelchairs, as well as the relevant normative references of this standard, in particular ISO 7176-14 and EN 1021-2 (or equivalent ISO 8191-2), and thus complies with EU Medical Device Regulation 2017/745 (MDR) and is classified as a Class I medical device under this regulation.

This product bears the UKCA marking in accordance with Part II UK MDR 2002 (as amended) Class I.


Medical device

For disposal of the Erivo and its components, see Chapter 12.5

Protect the device from moisture.

Indicates the temperature range within which the device can be used.

Observe the instructions for use!

These instructions for use contain, amongst other things, instructions, information and warnings regarding the operation of the Erivo and the charging of the batteries. These must be read and taken into account before putting the Erivo into service or charging it for the first time.

4


Indication of manufacturing date on the system label (see also Chapter 18).

Name and address of the device manufacturer (see reverse of this IFU).

Use the charger only indoors, not outdoors.

Double-insulated, equipment class 2

Risk of electric shock. Do not open the charger housing. Do not use the charger if cables are damaged.

Warning: risk of crushing hands, especially when folding and unfolding the Erivo — see Chapter 11.1.

Explanations of stickers on the Erivo: see Chapter 18.

1.7 Permitted Operating Conditions / Locations

- Observe the notes on permitted operating conditions in these instructions for use.
- Restrictions on permitted operating conditions (e.g. maximum gradient capability, maximum permitted obstacle height, maximum user weight, etc.) must be strictly observed!
- The Erivo may only be operated at temperatures between -10 °C and $+50\text{ °C}$. Do not expose the Erivo to heat sources (such as intense sunlight), as surfaces may heat up considerably.
- Avoid driving on unsurfaced ground (e.g. loose gravel, sand, mud, snow, ice or through deep puddles), as well as in rain, dense fog, thunderstorms, storms and poor visibility.
- Do not leave the Erivo in a damp or wet location for an extended period.
- In particular, observe the safety and hazard notices from Chapter 3 onwards.

The Erivo may rust, or electrical and mechanical defects may occur if rain, smog or moist air acts directly on the Erivo.

Do not expose the Erivo to prolonged intense sunlight when not in use. This causes the motor to heat up and in extreme cases it will not deliver its full power. Plastic parts also age more quickly under intense sunlight.

Performing so-called „wheelies“ (anti-tip wheels removed from the Erivo, drive wheels on the ground, front wheels (castors) hanging freely in the air) is not permitted. Alber GmbH disclaims all liability for accidents that occur as a result of this or similar driving behaviour.

1.8 Standard Scope of Delivery

- Erivo wheelchair
- Control unit
- Battery pack
- Charger
- These instructions for use

1.9 Optional Accessories

For information on available accessories, please refer to the Alber website.

It is recommended to use only original Alber accessories.

Alber is not liable for damage to the product or accidents (such as fires etc.) caused by malfunction of non-original accessories or spare parts. The warranty does not cover repairs required due to malfunction of non-original accessories. However, you may commission such repairs on a chargeable basis.

2. Commissioning

Your Erivo is delivered ready for operation with all standard operating components and any optionally ordered accessories.

Upon delivery of the Erivo, your specialist dealer will instruct you in the operation of the system and any accessories ordered. You will also be given these instructions for use, which contain, in addition to technical information, important notes on driving. Read the instructions for use carefully, in particular all safety and hazard notices.

The Erivo may only be driven for the first time after reading these instructions for use.





Before first use of the Erivo, its battery pack [7] must be fully charged (see Chapter 4.3).

In the event of a system failure or other serious system errors, the Erivo is brought to a halt quickly, as this normatively represents the safe state. The wheelchair user must be physically capable and have the reaction time required to counteract the resulting braking deceleration forces. If they are unable to do so, journeys with the Erivo must be carried out with the safety belt fastened.

Do not smoke while using the wheelchair. Burning cigarettes could fall unintentionally onto the seat cushion [14] and scorch or set it alight.

2.1 Adapting the Erivo to the User

2.1.1 Swapping the Armrests

For left-handed persons or physically impaired persons, it may be advantageous to swap the front parts [3a] of the armrests [3] - i.e. to have the control unit [2] fitted on the left instead of the right side, and the left front armrest part (with USB charging ports and lighting) fitted on the right side. Alber recommends having the swapping of the armrests [3] carried out by one of its representatives or an authorised specialist dealer.

Proceed as follows if you wish to swap the armrests:


- Switch off the electric wheelchair, see Chapter 5.2.
- Fold the armrests [3] upwards.
- Loosen three Allen bolts [58] on the underside of the rear part [3b] of each armrest [3] (see illustration) using an Allen key, key width 3.
- Remove the Allen bolts [58] and store them safely.
- Carefully pull the front parts [3a] of the armrests [3] out of the rear, fixed parts [3b].
- Disconnect the connection between the sockets [61] on the front parts [3a] and the plugs [60] on the rear parts [3b]:
 - Hold the front part of the first armrest [3a] firmly with one hand and simultaneously press the locking/unlocking tab [59] on the plug [60] of the rear part [3b].
 - Carefully pull the front part [3a] of the first armrest away from the fixed rear part [3b].
 - Repeat the above steps for the second armrest [3].
- Swap the front parts [3a] of the armrests.
- Connect the sockets [61] on the front parts [3a] to the plugs [60] on the rear parts [3b]:
 - Hold the front part [3a] of the first armrest firmly with one hand and carefully press the plug [60] into the socket [61] until you hear a quiet click.
 - Repeat this step for the second armrest [3].
- Carefully insert the front parts [3a] into the rear parts [3b].
Make sure that the cables are not trapped between the rear parts [3b] and front parts [3a]!
- Insert three Allen bolts [58] into each rear part [3b] and tighten them using an Allen key, key width 3 (tightening torque: 1.5 Nm).
- Fold the armrests [3] back down.

2.1.2 Extending / Shortening the Armrests

For persons with long arms, it may be more comfortable for the armrests to be longer than the standard setting. Alber recommends having the armrests extended by one of its representatives or an authorised specialist dealer.

Proceed as follows to extend the armrests [3]:

- Switch off the electric wheelchair, see Chapter 5.2.
- Fold the armrests [3] upwards.
- Loosen three Allen bolts [58] on the underside of the rear part [3b] of each armrest [3] using an Allen key, key width 3.
- Remove the Allen bolts [58].
- Carefully pull the front parts [3a] out of the rear, fixed parts [3b] far enough to allow two Allen bolts [58] to be inserted.
- Insert two Allen bolts [58] into each rear part [3b].
- Tighten with Allen key, key width 3 (tightening torque: 1.5 Nm).
- Store the removed Allen bolts [58] safely.

Proceed as follows to shorten the armrests [3]:

- Switch off the electric wheelchair, see Chapter 5.2.
- Fold the armrests [3] upwards.
- Loosen two Allen bolts [58] on the underside of the front part [3a] of each armrest [3] using an Allen key, key width 3, and remove them.
- Carefully push the front parts [3a] into the rear, fixed parts [3b] far enough to allow three Allen bolts [58] to be inserted.
- Insert three Allen bolts [58] into each rear part [3b] and tighten them with an Allen key, key width 3 (tightening torque: 1.5 Nm).
- Fold the armrests [3] back down.

2.2 Fitting / Adjusting the Safety Belt

The safety belt prevents the user of the Erivo from sliding out of the electric wheelchair when it comes to a standstill during travel.


The decision to fasten the safety belt is the responsibility of the user or the person responsible. For certain medical conditions, there is a medical indication to wear the safety belt.

When fitting the safety belt:

- Adjust the length of the safety belt to your body dimensions.

2.3 Attaching / Removing the Bag

Behind the backrest of the electric wheelchair there is a bag [6] in which you can store items needed during or after travel. The bag [6] has a capacity of 17 litres and a carrying capacity of 10 kg. The bag [6] is removable.



Proceed as follows to remove the bag [6]:

- Hold the bag [6] by its handle [6a].
- Pull the tab [6b] upwards, located behind the handle [6a] / below the push handle [5].
- The bag fastening releases and you can remove the bag [6].

Proceed as follows to attach the bag [6]:

- Hold the bag [6] by its handle [6a].
- Place the lower end of the bag [6] with its two lugs [6c] against the wheelchair (see illustration).
- Press the upper end of the bag [6] with your hand against the locking mechanism until it clicks audibly into place. The bag [6] is now secured to the wheelchair.


2.4 Using the Storage Basket

Beneath the seat cushion [14] there is a storage basket [11] in which you can place items needed during or after travel.

- When using the storage basket, note that it has a maximum load capacity of 5 kg.

2.5 Switching the System On


Before beginning a journey:

- Make sure that the Erivo is fully unfolded and that it can no longer be folded without operating lever [47] (see Chapters 11.1, 11.1.1 and 11.1.2).
- Make sure that the battery pack is fully charged (see Chapter 4.3).
- Make sure that you are correctly seated on the seat cushion [14].

If so, proceed as follows:

- Press the On/Off / Anti-theft lock button [20] on the control unit [2].
When the Erivo is switched on, the anti-theft lock is automatically activated and a padlock symbol [29] is visible on the screen [23] (see illustration).
After switching on, the red indicator on the battery pack [7] lights up.
- Press the On/Off / Anti-theft lock button [20] again to deactivate the anti-theft lock, i.e. to unlock the Erivo.
The padlock symbol [29] then disappears and the standard screen (see illustration and Chapter 5.13) is displayed.
- Begin your journey.
- If a journey cannot be started due to a fault or malfunction, this will be indicated on the screen [23] of the control unit [2] (see Chapter 6).

After switching off the system (see Chapter 5.2), both parking brakes [8] are automatically activated.

If you are not yet familiar with the Erivo and its driving characteristics, use Beginner Mode to first reduce the maximum speed to a minimum (see Chapter 5.8).

When driving the Erivo, observe the safety and hazard notices in Chapter 3 without fail.

The Erivo may only be operated at temperatures between $-10\text{ }^{\circ}\text{C}$ and $+50\text{ }^{\circ}\text{C}$. Do not expose the Erivo to heat sources (such as intense sunlight) as surfaces may heat up considerably.

Do not expose the Erivo to prolonged intense sunlight when not in use. This would cause the motor to heat up and in extreme cases it would not deliver its full power. Plastic parts also age more quickly under intense sunlight.

If you wish to push the Erivo manually, use Push Mode, see Chapter 5.10.

2.6 Switching the System Off

When you have finished your journey and no further journey will follow for an extended period, you should always switch off the Erivo. This saves energy and prevents the Erivo from being unintentionally set in motion by accidental contact with the joystick [24]. Please continue reading in Chapter 5.2.

2.7 Range

For every user of the Erivo, one of the most interesting questions is the system's range. In general, the range with the 12 Ah battery is approximately 13.4 km; with the 20 Ah battery approximately 21.5 km. These are ideal values and relate to flat, firm terrain. Deviations arise due to topographical conditions, ambient temperature, road surface, rear-wheel tyre pressure, frame geometry, wheelchair weight, user weight and tyre type.

3. Safety and Hazard Notices for Driving with the Erivo

3.1 General Notes

- The Erivo wheels are controlled via the joystick on the control unit (see Chapter 5.3). When the system is switched on, every joystick input is converted into a driving command, including forward/reverse movement, acceleration and braking. Switch off the control unit during any longer stop to prevent unintentional movement of the wheelchair caused by accidental contact with the joystick.
- It is recommended to carry out a short driving training session after receiving the Erivo, to familiarise yourself thoroughly with the drive and its capabilities.
- When transferring in or out of the Erivo, its drive wheels must be set to electric motor mode to prevent unintentional movement.

3.2 Driving Training Notes

- The safety and wellbeing of the driver come first. It is essential to become familiar with the driving characteristics of the Erivo. Your specialist dealer or Alber district managers will support you during a free device induction.
- Observe the information, safety and hazard notices in these instructions for use.
- Begin your first driving attempts with particular care on a flat, hard surface (e.g. asphalt, concrete, solid floor, not gravel etc.).
- Carry out your driving training in areas free of obstacles and pedestrians.
- Always adapt your speed to prevailing conditions in order to safely manoeuvre around suddenly appearing obstacles or bring the Erivo to a halt.
- When starting off, never push the joystick immediately to the full deflection. This may result in uncontrollable wheelchair movements (accident risk).
- Move the joystick smoothly and not jerkily. Avoid sudden full deflections, especially in potential hazardous situations you wish to avoid. Cautiously brake the Erivo to a standstill first.
- When you release the joystick, the wheelchair brakes gently. If emergency braking (immediate standstill) is required, perform an emergency stop (see Chapter 5.5).
- Brake the Erivo exclusively via joystick control. Never grab the tyres „for support“.
- Before driving the Erivo on descents or inclines, you should be confident in handling it on flat terrain.
- Always drive on inclines with a fully charged battery pack.
- Do not drive on descents with a fully charged battery pack, as this can cause damage to the battery pack and motors.
- Exercise special caution when driving on any type of descent.
- Drive on descents at a slow speed, especially to be able to swerve around sudden obstacles or stop the wheelchair in time.
- Never drive crosswise across steep slopes. An inadvertent shift in the centre of gravity could cause the wheelchair to tip sideways.
- When driving along sloping kerbstones etc., steer slightly into the slope. Cross smaller obstacles such as kerbstones never in parallel travel, but at right angles. This means both wheels cross the obstacle simultaneously, not staggered. Use only a low speed appropriate for the obstacle.
- Never remove the anti-tip wheels, not even to cross larger obstacles.
- When driving on public roads, observe the applicable road traffic regulations. Your wheelchair must be fitted with the additional equipment prescribed therein.

=> continue to the next page

**Caution when driving on descents with a fully charged battery pack!**

With a fully charged battery pack and high speeds, a reduction in speed or automatic shutdown of the system is possible. Therefore reduce your speed.

Overcome obstacles (e.g. kerbstones) as slowly and carefully as possible. The maximum permitted obstacle height is 4 cm. It is your discretion whether to seek assistance from another person.

In the event of a technical problem or error message, contact your specialist dealer immediately.

Observe the following safety and hazard warnings.

If possible, do not drive in reverse. Turn around instead and drive forward. If reversing cannot be avoided, drive slowly and with extreme caution. Additionally, use a helper who guides you and alerts pedestrians to your presence.

3.3 Safety Instructions

- The system must be switched off at the control unit before any work on the wheelchair begins.
- For safety reasons, avoid stopping on inclines and descents to prevent unintentional movement.
- In the (unlikely) event of overheating or fire of the battery pack, it must under no circumstances come into contact with water or other liquids. Battery manufacturers recommend sand as the only effective extinguishing agent.
- Using the Erivo as a vehicle seat is not permitted.

Before starting a journey:

- Driving the wheelchair with the anti-tip wheels removed is not permitted.
- Before every journey, check the condition of the Erivo wheels. If the tyres have reached their wear limit (tread depth less than 1 mm), the Erivo must no longer be operated.
- Regularly check the tyre pressure of the Erivo drive wheels. Correct pressure information is printed on the outside of the tyres and in the table in Chapter 17. Insufficient tyre pressure can strongly affect driving behaviour and range.
- Both Erivo drive wheels should always have the same tyre pressure.
- Check that fixings, connectors and removable parts are correctly fastened and connected (see Appendix B, section „Mechanics“).
- The Erivo must not be used in conjunction with a wheelchair towing device, handbike or other towing device.
- Never attach objects such as bags to the frame or armrests of the Erivo to avoid damage.
- The decision to fasten the safety belt is the responsibility of the user or the responsible person. For certain medical conditions, there is a medical indication to wear the safety belt.

During travel with the Erivo:

- Before driving the Erivo on descents or inclines, be confident in handling it on flat terrain.
- The maximum permitted gradient specified by the wheelchair manufacturer must not be exceeded.
- Test the brakes at the start of each journey. (Does the Erivo brake to a standstill when you release the joystick [24]?)
- Always begin journeys at the lowest speed level and increase it gradually according to external conditions.
- If possible, do not drive in reverse. Turn around and drive forward instead. If reversing cannot be avoided, drive slowly and with extreme caution. Use a helper to guide you and alert pedestrians.
- Approach stairs or drop-offs with extreme caution.
- When the Erivo is switched on, every touch of the joystick is converted into a driving command. Remove your hand from the control unit when stopping or waiting at potential hazardous locations (e.g. at a pedestrian traffic light, on inclines and descents, or at ramps of any kind).
- Never reach into the wheel while the Erivo is switched on.
- Do not insert or throw objects of any kind into a switched-off or switched-on Erivo wheel.
- Never attach objects of any kind to the Erivo wheels! This could cause damage.
- When driving on pavements, maintain sufficient clearance (preferably at least one wheelchair width) from the kerb edge.
- Avoid driving on unsurfaced ground (e.g. loose gravel, sand, mud, snow, ice or through deep puddles).
- Avoid driving over sharp or pointed objects such as sharp edges, stones or broken glass.
- Never leave the Erivo unattended, whether switched on or off.

- The Erivo can be affected by strong electromagnetic fields.
- In rare cases, the Erivo may interfere with other installations, e.g. anti-theft systems in department stores.
- Escalators and travelators must not be used with the Erivo.
- Never suddenly change direction left or right while driving at maximum speed, as this may cause the wheelchair to tip sideways.
- Reduce your speed when cornering and when driving in confined spaces.
- If you intend to stop on an incline or descent, the wheelchair must be positioned crossways to the slope.
- Never drive crosswise on descents.
- Never reach into the Erivo drive wheels or front wheels while driving.
- Brake the Erivo smoothly and in proportion to your speed (not jerkily).
- When travelling in vehicles, sit exclusively on the installed seats with the corresponding restraint systems. Failure to do so may result in injury to you and other passengers in an accident.
- When travelling in vehicles, secure your Erivo in accordance with the applicable statutory provisions and regulations.
- Using the Erivo as a vehicle seat is not permitted.
- When driving on public roads and pavements in Germany, the Road Traffic Regulations (StVO) and the Road Vehicle Licensing Regulations (StVZO) must be observed. In other countries, the applicable national regulations must be observed.
- The Erivo is intended only for the transport of persons with restricted mobility and must not be misused, e.g. for children playing, carrying children, or transporting goods.
- Operation of the Erivo is not permitted in the vicinity of strong magnetic fields, such as those generated by permanent magnets, transformers, tomography equipment, etc.
- Avoid driving with the Erivo in adverse conditions, e.g. in storms, hail, and in dense undergrowth.

After travel with the Erivo:

- Switch off the Erivo immediately when not in use to prevent accidental triggering of a driving impulse through contact with the joystick, and to prevent self-discharge of the battery pack.
- Charge the battery pack of your Erivo after every journey where possible.

3.4 Special Safety Instructions for Slopes, Descents and Ramps

To achieve maximum wheelchair stability and thus maximum driver safety, we recommend the following seating positions:

Travel on flat ground and on descents

Travel on inclines and ramps

Maximum gradient/descent: 9 degrees
Maximum user weight: 150 kg

The following points must also be observed without fail:

- The maximum user weight is 150 kg.
- Only gradients and descents of maximum 9 degrees may be driven. Never drive gradients exceeding the rated gradient.
- Drive on inclines, descents and in particular ramps exclusively at a low speed appropriate for the situation.
- Avoid tight turns on inclines and descents.
- In particular on ramps, do not drive curves.
- Avoid driving on long ramps and in particular on uneven terrain, to avoid the risk of motor failure.
- The ability of your Erivo to negotiate inclines is influenced by, among other things, your weight, speed and approach angle.
- Observe the seating positions illustrated above, especially on inclines and ramps.

3.5 Obstacles

- The maximum permitted obstacle height is 4 cm.
- Caution! Drive slowly and carefully towards an obstacle until the Erivo wheels touch it. Then carefully overcome the obstacle. It is your discretion whether to seek assistance from another person.

3.6 Hazardous Locations and Situations

- The wheelchair user decides independently and on their own responsibility, considering their driving skills and physical capabilities, which routes to drive.
- Before starting a journey, they must check the Erivo wheels for worn or damaged tyres, the battery pack charge level, and the functionality of the warning signal (horn) on the control unit.
- These safety checks and personal driving skills are of particular importance at the following hazardous locations, where driving is at the discretion and own risk of the Erivo driver:
 - Quay walls, landing and mooring stages, paths and places near bodies of water, unsecured bridges and dikes
 - Narrow paths, descents (e.g. ramps and driveways), narrow paths on slopes, mountain routes
 - Narrow and/or sloping paths along main and secondary roads or near drop-offs
 - Routes covered with leaves, snow or ice
 - Ramps and lifting devices on vehicles
 - Cobbled paths/roads

When cornering or turning on uphill or downhill gradients, a shift in the centre of gravity may increase the risk of the vehicle tipping over. You should therefore carry out these manoeuvres with extra care and only at low speed!

When crossing roads, intersections and level crossings, heightened caution is required. Never cross rails in roads or at level crossings in parallel travel, as the castors of the wheelchair could become jammed.

Exercise particular caution when using ramps and lifting devices on vehicles. During the lifting or lowering of the ramp or lifting device, the Erivo must be parked and switched off with sufficient clearance from the edges.

Do not drive directly up to or to the edge of a ramp or lifting device, maintain sufficient clearance.

In wet conditions, tyre grip on the surface is reduced and there is an increased risk of slipping. Adapt your driving behaviour accordingly and never drive with worn tyres.

The Erivo should not be used in rain or snowfall to avoid technical damage.

Be aware that when driving on descents, the braking distance of the Erivo can be significantly longer than on flat terrain, depending on speed and gradient. Adjust your speed accordingly.

Do not drive on descents with a fully charged battery pack, as this can cause damage to the battery pack and motors.

Exercise heightened caution when driving on cobbled paths/roads, as the castors of the wheelchair may become jammed in the joint gaps.

In the event of fire or smoke, wheelchair users are particularly at risk of being killed or seriously injured as they are unable to move away from the source of fire or smoke. Lit matches, lighters and cigarettes can cause an open flame in the vicinity of the wheelchair or on clothing.

Avoid using or storing the wheelchair near open flames or flammable products, and do not smoke while using the wheelchair.

Using the Erivo as a vehicle seat is not permitted.

Performing so-called „wheelies“ (anti-tip wheels removed, drive wheels on the ground, front castors hanging freely in the air) is not permitted. Alber GmbH disclaims all liability for accidents resulting from this driving behaviour.

3.7 Electromagnetic Interference (EMI)

Any device with electrical components can cause interference and affect the operational capability of other devices. Your Erivo may therefore be affected by radio and television transmitters, amateur radio transmitters, radios and mobile phones. In extreme cases, this could cause the Erivo to release its brakes and move on its own. The Erivo's control system may also be permanently damaged.

Your Erivo can withstand electromagnetic interference of up to 20 V/m.

EMI Warnings

- In general, travel near electrical interference fields should be avoided.
- Do not use portable transceivers such as walkie-talkies or mobile phones while the Erivo is switched on.
It is better to switch these devices off completely before a journey and switch them back on only after the journey has ended.
The Erivo should be completely switched off before operating such devices.
- Be aware of nearby powerful transmitters such as radio or television stations and avoid getting close to them.
- Switch off the Erivo completely if you notice unintentional, independent movement of the Erivo during electric motor operation. Re-commission it only in an environment where there is no or only minimal electromagnetic interference.
- Note that adding accessories and components, as well as any modification, may make the Erivo more susceptible to electromagnetic interference.

In rare cases, operation of the Erivo may interfere with other installations, e.g. anti-theft systems in department stores.






4. Battery Pack

Before commissioning the Erivo and before starting the charging process, read and observe the general information and notices, as well as the safety and warning notices in Chapters 4.5 to 4.10. Failure to comply with the safety notices and instructions may damage the product or result in electric shock, fire and/or serious injury.

Alber GmbH accepts no liability for damages resulting from failure to observe the general notices and safety instructions.

Use only the battery pack supplied with the Erivo or one supplied by Alber. Use of third-party products will void warranty and guarantee entitlements. Alber also accepts no liability for any resulting damage.

4.1 Removing the Battery Pack

As a rule, the battery pack [7] will always remain on your wheelchair. If it needs to be removed (e.g. for transport), proceed as follows:

- Switch off the system at the control unit (see Chapter 5.2).
- Turn the knobs [45] on both sides of the battery pack [7] from the horizontal (locked) position 90 degrees to the vertical (unlocked) position (see the illustration and the text printed on the battery pack [7]).
When the knobs [45] are turned 90 degrees, you will hear a click.
- Turn the battery clamps [46] on both sides of the battery pack approximately 90 degrees away from the battery pack [7] (see illustration).
- Remove the battery pack [7].
- Store the battery pack [7] in a safe place.

4.2 Fitting the Battery Pack

If the battery pack has been removed from the wheelchair, it must be refitted to operate the Erivo in electric motor mode.

- Make sure that there are no foreign particles (especially metallic objects) in the system power socket [61] or in the charging socket [49] on the battery pack. Clean these if necessary.
- Make sure the left and right battery clamps [46] are approximately transverse to the direction of travel.
- Place the battery pack [7] into the battery compartment so that it rests against the vertical (rear) part.
- Turn the battery clamps [46] on both sides towards the battery until you hear a click (see illustration).
The clamps [46] are now locked.
- Turn the knobs [45] on both sides from the vertical (unlocked) to the horizontal (locked) position (see the illustration and the text printed on the battery pack [7]).
When the knobs are turned 90 degrees, you will hear a click.

4.3 Charging the Battery Pack

After each extended use of your Erivo, its battery pack should be recharged. This can be done via the control unit (Chapter 4.3.1) or directly at the removed battery pack (Chapter 4.3.2).

The charging process may take up to approximately 4.5 hours for the 12 Ah battery and up to 7.5 hours for the 20 Ah battery.

Do not interrupt the charging process before the battery pack [7] is fully charged. The LED indicator [51] on the charger will then light up green.

The battery pack [7] should always be fully charged!

The battery pack [7] may only be charged indoors and never outdoors. Avoid direct sunlight during charging (e.g. through a window) to prevent the charger from overheating. Before charging the battery pack [7], observe the following:



Only an Alber-supplied charger [50] may be used to charge the battery pack [7]. Using a different charger is strictly prohibited!

The charger [50] must not remain permanently connected to the Erivo. Disconnect it from the mains socket and then from the battery pack [7] as soon as charging is complete.

Always charge the battery pack [7] fully. Only interrupt the charging process in an emergency. The duration of charging depends on the remaining capacity of the battery pack [7].

4.3.1 Charging at the Control Unit

You can charge the battery pack [7] whether the Erivo is switched on or off.

To charge the battery pack [7] on the control unit [2], proceed as follows:

- Switch off the Erivo if desired (see Chapter 5.2).
- Insert the charging plug of the charger [50] into the charging socket [49] under the control unit [2].
- Insert the mains plug [53] of the charger [50] into a suitable socket. The LED indicator [51] on the charger [50] will light up red. This indicates that charging has started but is not yet complete.


If the Erivo is switched on, „charging“ [40] will be displayed on the screen [23] during charging, and the headlight on the battery pack [7] will flash alternately brighter and dimmer.

- When charging is complete, the LED indicator [51] on the charger [50] lights up green.


If the Erivo is switched on, „full charged“ [41] will appear on the screen after charging is complete.

After charging, disconnect the mains plug [53] from the socket.


- Disconnect the charging plug [52] from the charging socket [49] under the control unit.

During the charging process, the control unit [2] is deactivated; electric motor operation and joystick control are not possible.


The parking brakes [8] are automatically activated when the Erivo is switched on and the battery pack [7] is being charged at the control unit [2].

4.3.2 Charging at the Battery Pack Charging Socket

Follow these steps to charge the battery pack once you have removed it from the wheelchair:

- Switch off the Erivo (see Chapter 5.2).
- Remove the battery pack [7] from the wheelchair (see Chapter 4.1).
- Insert the charging plug [52] of the charger [50] into the charging socket [49] of the battery.
- Insert the mains plug [53] into a suitable socket.
The LED indicator [51] on the charger [50] lights up red. This indicates that charging has started but is not yet complete.
When charging is complete, the LED indicator [51] on the charger lights up green. After charging, disconnect the mains plug [53] from the socket.
- Disconnect the charging plug [52] from the charging socket [49].
- Refit the battery to the electric wheelchair (see Chapter 4.2).

4.4 Capacity Indicator

The battery pack [7] required to operate the Erivo is fitted as standard behind the seat and has no charge indicator of its own. The available battery capacity is therefore shown on the screen [23] of the control unit. This display includes a percentage indicator [30] and a graphical bar indicator within a stylised battery icon [31].

4.5 General Charging Notes

- Never charge the battery pack [7] in the presence of or near flammable liquids or gases.
- In the (unlikely) event of overheating or fire of the battery pack [7], it must under no circumstances come into contact with water or other liquids. Battery manufacturers recommend sand as the only effective extinguishing agent.
- Never carry out the charging process in rooms where moisture could condense on the battery pack [7].
- Carry out charging with the wheelchair in a location that offers at least twice the floor space of the wheelchair and adequate ventilation so that there is no risk from accumulation of flammable gas.
- During charging, explosive gases may be produced. Keep the wheelchair and charger away from ignition sources such as flames and sparks.
- Never charge the battery pack [7] at temperatures below 0 °C or above 40 °C.
- Under no circumstances use an extension cable to connect the charging plug [52] to the charging socket [49] of the battery pack [7] or the charging socket [49] of the control unit [2].

4.6 Safety Instructions for the Charger and Charging Process

- When the battery pack [7] is charged directly on the wheelchair, it must be secured against unintentional rolling by activating electric motor mode (see Chapter 5.1)!
- Use only the supplied charger [50] to charge the battery pack [7]. The charging process ends automatically when the battery pack [7] is charged, so overcharging is excluded.
- Never use a charger not supplied by Alber.
- Do not charge any other batteries with the charger; only charge the Erivo's battery pack [7].

- The charger [50] must not remain permanently connected to the Erivo. Check charging progress periodically and disconnect the charger from the mains socket, then from the battery pack [7], as soon as charging is complete.
- Do not expose the charger [50] to any moisture (water, rain, snow) during the charging process.
- Caution regarding condensation. If the charger [50] is brought from a cold into a warm room, condensation may form. In this case, delay use of the charger [50] until the condensation has evaporated. This may take several hours.
- Never carry out the charging process in rooms where moisture could condense on the charger [50].
- Never carry the charger [50] by the mains cable or charging cables.
- Never pull the mains cable to disconnect the charger [50] from the socket.
- Route the mains cable and charging cable so that no one can step on or trip over them. Do not subject either cable to other harmful influences or loads.
- Do not operate the charger [50] if the mains cable, charging cable or plugs [52] [53] are damaged. Damaged parts must be replaced immediately by an Alber-authorized specialist dealer.
- Never use or dismantle the charger [50] if it has suffered a heavy blow, been dropped, or otherwise damaged. Bring the damaged charger [50] to an Alber-authorized medical supply dealer for repair.
- The charger [50] must not be used by young children.
- The charger [50] may only be operated with 100 V – 240 V AC mains voltage.
- Do not dismantle or modify the charger [50].
- Do not cover the charger [50] during charging and do not place objects on the device.
- Do not subject cables and plugs [52] [53] to pressure. Strong pulling or kinking of cables, trapping cables between a wall and a window frame, or placing heavy objects on cables or plugs could lead to electric shock or fire.
- Never short-circuit the pins of the charging cable plug [52] with metallic objects.
- Make sure the mains plug [53] is firmly inserted in the socket.
- Do not touch the mains plug [53] or charging cable plug [52] with wet hands.
- Do not use the charger plug [52] and/or mains plug [53] if they are wet or dirty. Clean them with a dry cloth before use.
- After charging, first disconnect the charger mains plug [53] from the socket, then the charging plug [52] from the socket on the battery pack [7].
- Ensure that no moisture can enter the socket [49] on the battery pack [7] after disconnecting from the charger [50].
- Store the charger [50] in a dry, well-ventilated environment.
- The charger [50] must not be stored in locations where it is exposed to direct heat sources, moisture, water, oil, acids or alkalis, or where ozone may be generated.

The Alber-supplied chargers [50] may only be used to charge the Erivo's battery pack [7]. These chargers must not under any circumstances be used to charge other batteries.

Before and after the charging process, ensure that the socket [49] on the battery pack is clean and free of metallic particles. If any are present, they must be removed with a clean, dry cloth!

4.7 Safety Instructions for the Battery Pack

- Before first use of the Erivo, its battery pack [7] should be fully charged.
- The Erivo may only be operated at temperatures between $-10\text{ }^{\circ}\text{C}$ and $+50\text{ }^{\circ}\text{C}$.
- Never charge the battery pack [7] at temperatures below $0\text{ }^{\circ}\text{C}$ or above $45\text{ }^{\circ}\text{C}$.
- The battery pack [7] must not be exposed to heat (e.g. radiators) or fire. External heat exposure can cause the battery pack [7] to explode.
- In the (unlikely) event of overheating or fire of the battery pack [7], it must under no circumstances come into contact with water or other liquids. Cell manufacturers recommend sand as the only effective extinguishing agent.
- Your Erivo consumes energy with every use. Therefore, charge the battery pack [7] fully after each use where possible.
- The housing of the battery pack [7] must not be opened or dismantled. Improper opening or deliberate destruction poses a risk of serious injury. Opening also voids any warranty claims.
- Never connect the socket [61] of the battery pack [7] to metallic objects, and ensure that the contacts do not come into contact with metallic objects (e.g. metal swarf) under any circumstances.
- If the socket [61] is dirty, clean it with a clean, dry cloth.

=> continue to the next page

- Never immerse the battery pack [7] in water.
- If the battery pack [7] is damaged or defective, it must be checked. Please contact your specialist dealer and clarify the procedure for return and repair (see also Chapter 4.9).
- If the battery pack [7] is defective or damaged, the Erivo must not be used further. Contact your specialist dealer regarding repair.
- Always keep the battery pack [7] clean and dry.

4.8 Storing the Battery Pack

- The lifespan of a battery pack [7] depends in part on how it is stored. Do not leave the battery pack [7] and Erivo wheels in hot locations for an extended period. In particular, car boots of vehicles parked in the sun should only be used for transport, not as a general storage location.
- Store the battery pack [7] in a cool, dry and well-ventilated environment, protected from damage and unauthorised access.
- To achieve optimum battery pack lifespan, it should be stored at a temperature of n.a. to n.a., an air pressure of n.a. to n.a. hPa and a humidity of n.a. to n.a. %. Under these conditions the storage life of the battery pack [7] is 6 months.
- Before storage, the battery pack [7] must be fully charged.
- If the battery pack [7] is stored for longer than 6 months, it must also be fully charged every 6 months.
- Do not expose the battery pack [7] to any moisture (water, rain, snow, etc.) during storage.
- Charge the battery pack [7] before storing and check the charge level every month. Recharge if the charge level falls below 50%.

4.9 Safety and Warning Notices for Transport and Shipping

- The Erivo battery pack [7] uses lithium-ion cells. Corresponding statutory regulations therefore apply to transport and shipping of the battery pack [7] and must be strictly observed! For example, a defective battery pack [7] must under no circumstances be transported in aircraft, either as hand luggage or as checked luggage.
- If your battery pack [7] is defective, please bring it personally to your specialist dealer, as postal or courier shipping is also strictly regulated for lithium-ion batteries. We also recommend contacting the specialist dealer in advance.
- Since transport regulations may change annually, we strongly recommend contacting your tour operator or airline/shipping company before travel and obtaining information about the currently applicable regulations.

If you require packaging compliant with statutory regulations for transport purposes, you can request this from your specialist dealer.

Discuss the necessary measures with your specialist dealer before shipping the battery pack [7].


4.10. Important Information for Air Travel

The Erivo, and in particular the lithium-ion battery contained within it, is technically designed in accordance with applicable standards, directives and regulations. It meets the requirements of the *International Air Transport Association (IATA)* for carriage by air. We confirm this annually with our transport certificate, which is prepared on the basis of each year's new edition of the IATA Dangerous Goods Regulations (DGR), which typically appears in mid-December. The current certificate is available for download from the Alber website (<https://www.alber.de>); we are also happy to send it to you on request.

Nevertheless, when planning a trip, contact your tour operator and the airlines regarding taking the Erivo on all flights booked. Ultimately, it is at the pilot's discretion whether to accept the Erivo on board or to decline its transport. Neither Alber GmbH nor its medical supply dealer can influence this decision.

When travelling, also observe the applicable regulations for the transport of lithium-ion batteries in the countries you visit.

5. Driving the Erivo

5.1 Switching the System On

Do not touch the joystick [24] during the switch-on process. If you do, an error message will be generated.

- Briefly press the On/Off / Anti-theft lock button [20]. The system switches on. The start-up screen is first displayed (see Chapter 5.13).
- Press the On/Off / Anti-theft lock button [20] again to cancel the anti-theft lock and drive the Erivo (see Chapter 5.12). An open padlock symbol briefly appears, indicating the anti-theft lock is cancelled.

The standard screen (see Chapter 5.13) is then displayed, and the headlight on the battery pack [7] lights up red.

- If you are not yet familiar with the Erivo and its driving characteristics, use Beginner Mode (see Chapter 5.8) or the speed level buttons [26] [27] to first reduce the maximum speed to a minimum (speed levels 1 and 2) (see Chapter 5.7).
- You can now drive the Erivo using the joystick [24]. However, first read the safety and hazard notices in Chapter 3.
- If the Erivo is not operated, it will automatically switch off after 30 minutes (see Chapter 7.2).


If a journey is not possible after switching on, there may be a fault. Chapter 6 contains a list of error codes.

5.2 Switching the System Off

When a journey is finished and you wish to leave the seat or no further journey will follow for an extended period, always switch off the Erivo. This saves energy and prevents the Erivo from being unintentionally set in motion by accidental joystick contact or by leaving the seat.

- Press and hold the On/Off / Anti-theft lock button [20] for approximately three seconds. The Erivo switches off and the screen [23] and the indicator on the battery [7] go out.

After switching off the system, both parking brakes are automatically activated.

5.3 Driving with the Joystick

The Erivo is driven with the joystick [24], via which the driver controls both speed and direction.

First, the maximum speed to be reached is preset using the speed level buttons [27] and [26]:

- Press the „higher speed level“ button [27] to increase the maximum speed.
- or:
- Press the „lower speed level“ button [26] to reduce the maximum speed.

The selected speed level is shown on the control unit screen as a black number on a white background (see Chapter 5.13).

For information on speed levels, see Chapter 5.7.

Begin driving as follows:

- Push the joystick [24] forward to drive forward.
- or:
- Push the joystick [24] backwards to drive backwards.

If possible, do not drive in reverse. Turn around and drive forward instead. If reversing cannot be avoided, drive slowly and with extreme caution. Use a helper to guide you and alert pedestrians.

=> continue to the next page

Accelerate as follows:

- Push the joystick [24] further forward or backward.

Brake gently as follows:

- Push the joystick [24] less far forward or backward.

Brake firmly as follows:

- Release the joystick [24] during travel.

The electric wheelchair will come to a standstill in the shortest possible time.

Steer as follows:

- Push the joystick [24] to the right to turn right.

or:

- Push the joystick [24] to the left to turn left.

Maximum speed is reached at full joystick deflection. During reverse travel, speed is automatically reduced for safety reasons.

Start your first driving attempts at a low speed to get used to the Erivo's handling characteristics.


Adapt speed to each driving situation. Reduce maximum speed especially in difficult driving conditions and in indoor areas to prevent accidents.

Hold on firmly to the armrests when the Erivo brakes.

Do not lean on the control unit — it could be damaged. Always use the wheelchair armrest to push yourself up, to stand, or to reposition yourself on the seat!

5.4 Using the Horn

The horn can be used to warn persons of hazards.

- Press the horn button [22] on the control unit [2].
The horn sounds.

5.5 Emergency Braking

In critical situations or in the event of unexpected Erivo movements, the emergency braking function can be used to bring the electric wheelchair to a very rapid halt.

- Press the On/Off / Anti-theft lock button [20] during travel to perform emergency braking.

During emergency braking, hold on firmly to the armrests [3] to prevent sliding out of the wheelchair.

Do not use the emergency braking function under normal conditions: frequent use of the emergency braking function can cause motor damage!

5.6 Triggering an Emergency Call

In an emergency (e.g. if you suddenly feel unwell), you can use the electric wheelchair's emergency call to alert other persons to your situation and call for help.

You can use the emergency call at any time when the Erivo is switched on, regardless of the driving mode or speed level, and even if the anti-theft lock is still active directly after switching on (see Chapter 5.1).

In an emergency, trigger an emergency call as follows:

- With the electric wheelchair switched on, hold down the emergency call button [21] for three seconds.



A voice announcement „Alarm“ then sounds for approximately five seconds and the SOS symbol [38] is displayed alternately white and red on the screen.

Use the emergency call exclusively in an emergency — never as a joke or in non-emergency situations.

Since the volume of the emergency call is relatively low, it can only alert persons in your immediate vicinity. In loud ambient conditions, the emergency call may be overlooked.

5.7 Limiting Maximum Speed (Speed Levels)

The minimum maximum speed of your Erivo is 1.3 km/h, which can be increased in individual levels up to a maximum speed of 6 km/h at full joystick deflection.

For various driving situations, for example in indoor areas or if you are not yet experienced with the Erivo, the maximum speed can be individually limited. Proceed as follows:

- Press the „lower speed level“ button [26].

or:

- Press the „higher speed level“ button [27].

The selected speed level is displayed on the control unit screen as a black number on a white background (see Chapter 5.13).

You can also use Beginner Mode to limit the maximum speed (see Chapter 5.8).

=> continue to the next page

Speed levels — Forward travel (see illustration, item number [33]):

- Speed level 1: You can drive up to 1.3 km/h.
- Speed level 2: You can drive up to 2.7 km/h.
- Speed level 3: You can drive up to 3.9 km/h.
- Speed level 4: You can drive up to 5.3 km/h.
- Speed level 5: You can drive up to 6.0 km/h.

The maximum speed set in this way is reached when full joystick deflection is reached.

During reverse travel, speed is automatically reduced for safety reasons.

Speed levels — Reverse travel (see illustration, item number [33]):

- Speed level 1: You can drive up to 0.9 km/h.
- Speed level 2: You can drive up to 1.1 km/h.
- Speed level 3: You can drive up to 1.4 km/h.
- Speed level 4: You can drive up to 1.6 km/h.
- Speed level 5: You can drive up to 1.7 km/h.

Speed levels can be set both while stationary and while driving.


Start your first driving attempts at a low speed to get used to the Erivo's handling.

Adapt speed to each driving situation. Reduce maximum speed especially in difficult driving conditions and indoors to prevent accidents.

If possible, do not drive in reverse. Turn around and drive forward. If reversing cannot be avoided, drive slowly and with extreme caution. Use a helper to guide you and alert pedestrians.

5.8 Driving in Beginner Mode / Exiting Beginner Mode

For various driving situations — e.g. if you are not yet familiar with the Erivo, or for driving indoors — the maximum speed can be limited. Select Beginner Mode for this:

- Switch on the Erivo (see Chapter 5.1).
 - Press the On/Off / Anti-theft lock button [20] again to cancel the anti-theft lock. The lock icon will then disappear and the default screen will be displayed (see the illustration, e.g. in Chapter 5.13).
 - Press the drive programme selector button [25]. The „Novice“ symbol [36] and the numbers 1 and 2 for speed levels 1 and 2 appear in the lower right of the screen. The number 1 is highlighted with a white background to indicate that speed level 1 is selected.
 - Push the joystick [24] forward to drive forward.
- or:
- Push the joystick [24] backward to drive backward.

If possible, do not drive in reverse. If unavoidable, drive slowly and with extreme caution. Use a helper to guide you and alert pedestrians.

If you wish to drive faster:

- Push the joystick [24] further forward.

and/or:

- Press the „higher speed level“ button [27] to select speed level 2.

If you wish to drive slower:

- Push the joystick [24] less far forward.

and/or:

- Press the „lower speed level“ button [26] to select speed level 1.

To exit Beginner Mode and return to general driving mode:



- Press the drive programme selector button [25]. The „Novice“ symbol [36] disappears. The standard screen with speed level display [33] is now visible (see Chapter 5.13).

Information:

- In Beginner Mode, speed levels 1 and 2 are available.
- For the maximum speeds achievable in forward and reverse travel at speed levels 1 and 2: see Chapter 5.7
- The headlight [1] in the armrest [3] automatically lights up while driving.

For braking and steering, observe the notes in Chapters 5.3 and 5.5.

5.9 Driving with Speed Control

If you do not wish to control the speed with the joystick [24], you can use Speed Control. The electric wheelchair will then drive at a fixed speed depending on the speed level:

- Speed level 1: 1.8 km/h
- Speed level 2: 2.6 km/h
- Speed level 3: 3.9 km/h
- Speed level 4: 5.3 km/h
- Speed level 5: 6.0 km/h

You can switch from a lower speed level to a higher one, or vice versa, whilst driving. You can also deactivate speed control whilst driving.

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- Press „lower speed level“ button [26] once or multiple times.

To end Speed Control:

- Push the joystick [24] backward
The speed control symbol [32] disappears and the Erivo brakes to a standstill.

or:

- Press the On/Off / Anti-theft lock button [20].
The padlock symbol [29] is displayed and the Erivo brakes to a standstill.

If, after deactivating speed control, you wish to continue driving using the on/off/immobiliser button [20] (electrically operated), proceed as follows:

- Press the on/off/immobiliser button [20] again.

When driving with speed control, the light [1] in the armrest [3] comes on automatically.

Before activating Speed Control, preselect speed level 1 so that the Erivo starts slowly. If it starts too fast, there is a risk of property damage and/or personal injury.

When driving with Speed Control, the joystick [24] no longer controls speed. You can only steer or end Speed Control with the joystick [24].

If Speed Control is activated but not started within five seconds, it is automatically cancelled and the speed control symbol [32] disappears.

5.10 Pushing the Erivo (Push Mode)

In addition to electric motor operation, the Erivo also offers the option of manual operation in Push Mode by an accompanying person. In this mode, the Erivo rolls particularly easily.

Activate Push Mode as follows:

- Switch on the system at the control unit [2] (see Chapter 5.1).
- Simultaneously hold down „higher speed level“ button [27] and „lower speed level“ button [26] for three seconds. This decouples the wheels and sets them to Push Mode. The Push Mode symbol [39] is then displayed on the screen [23].

In Push Mode, both drive wheels are in freewheel. This means no braking function is available. For longer stops, the wheels should therefore be switched back to electric motor mode to prevent unintentional movement of the Erivo.

Decoupling from electric motor mode to Push Mode must for safety reasons only be performed on flat ground and when stationary. The drive wheels should be unloaded when decoupling. Operating the clutch under load or during travel can damage the internal mechanism and is therefore not covered by warranty or guarantee. Decoupling on a descent or ramp can lead to a dangerous situation, as the Erivo may immediately or subsequently switch to freewheel and in the worst case move uncontrollably — it may rotate and/or roll downhill. Persons may be injured and/or objects damaged. Decoupling on a slope should only be done in emergency situations and only when an accompanying person is present who is able to hold the wheelchair back manually.

End Push Mode as follows:




- Move the joystick [24] in Push Mode. The Push Mode symbol [39] disappears, the anti-theft lock is automatically activated, and the start-up screen with the anti-theft lock symbol [29] is displayed.


To resume electric motor travel after exiting Push Mode:

- Press the On/Off / Anti-theft lock button [20]. The start-up screen will then appear (see Chapter 5.13).
- Press the On/Off / Anti-theft lock button [20] again to cancel the anti-theft lock (see Chapter 5.1)

To reactivate Push Mode after exiting it (with the Erivo still switched on):

- Simultaneously hold down „higher speed level“ button [27] and „lower speed level“ button [26] for three seconds. Dadurch werden die Räder ausgekuppelt und in den Schiebetrieb versetzt. The Push Mode symbol [39] is displayed on the screen [23] (see Chapter 5.13).


5.11 Releasing / Activating the Parking Brakes

To release one or both parking brakes [8] (e.g. when errors 40 to 42 are shown on the Erivo screen — see Chapters 5.12 and 6)):

- Use your foot to press down the red lever of the relevant parking brake [8] or the red levers of both parking brakes [8].

To activate one or both parking brakes [8]:

- Use your foot to push up the red lever of the relevant parking brake [8] or the red levers of both parking brakes [8].


5.12 Activating / Deactivating the Anti-Theft Lock

You can activate the anti-theft lock to prevent unauthorised persons from using the Erivo. To activate the anti-theft lock:

- Press the On/Off / Anti-theft lock button [20] when the Erivo is switched on. The screen shows the anti-theft lock padlock symbol [29] and the electric wheelchair cannot be operated.

To deactivate the anti-theft lock:

- With the anti-theft lock active, press the On/Off / Anti-theft lock button [20].. The padlock symbol [29] disappears and the electric wheelchair can be operated.


29

5.13 Screen Displays

Start-up Screen

After switching on the Erivo, the start-up screen with the anti-theft lock symbol [29] is displayed.

Standard Screen

The following displays are shown on the standard screen:

- Battery charge level display as percentage [30] (see also Chapter 4.4).
- Graphical battery charge display [31] (see also Chapter 4.4).
- Speed control display [32] (see also Chapter 5.9)
Visible when Speed Control is activated.
- Speed level display [33]
Selectable levels are shown in grey on black.
The selected level is shown in black on white.
- Bluetooth® connection indicator [34]
Visible when Erivo has an active Bluetooth® connection.
- Erivo speed display in km/h [35]

Beginner Mode

- When Beginner Mode is selected, the „Novice“ display [36] is shown (see Chapter 5.8).

Error Code Display


- Errors are displayed on the screen with the word „Error“ and a number or alphanumeric combination. For the meaning of error codes [37] and possible remedies: see Chapter 6. The illustration shown here is just an example of an error code!

Emergency Call Display


- The „SOS“ display [38] flashes when you have triggered an emergency call (see Chapter 5.6).


Push Mode Display

- The Push Mode display [39] is shown when Push Mode is selected (see Chapter 5.10).

Battery Charging Display

- The battery charging display [40] is shown when the Erivo is switched on while the battery pack [7] is being charged (see Chapter 4.3.1).

Battery Full Display


- The battery full display [41] is shown when the Erivo is switched on and the battery pack [7] is fully charged at the end of the charging process (see Chapter 4.3.1).

For information on Erivo malfunctions and possible remedies, see Chapter 6.

5.14 Using the Lighting

For safety reasons, it is recommended to switch on the lighting at dusk and in the dark.


To switch on the lighting:

- Touch the slightly recessed surface [1a] on the underside of the armrest for approximately one to two seconds. The headlight [1] illuminates.

To switch off the lighting:

- Touch the slightly recessed surface [1a] on the underside of the armrest for approximately one to two seconds. The headlight [1] goes out.

5.15 Using the USB Charging Ports

You can use the USB charging ports to charge an electrical device such as your smartphone or tablet. There is one USB-A port and one USB-C port in the armrest. Follow these steps if you wish to use one of the USB charging ports [56] on the electric wheelchair:

- Press the cover [55] of the USB ports [56] backward with one finger and hold it there.
- With the other hand, insert a charging plug into the appropriate USB port [56].
- Release the cover [55].
- After charging, remove the charging plug from the USB port [56].
- Close the cover [55] by pressing it downwards.

6. Faults

If errors other than those described in Chapter 6.1 are displayed on the screen [23], or other faults occur than those listed in Chapter 6.2, there is a serious system malfunction. The Erivo is then no longer operational. It must be repaired by your medical supply dealer or sent to Alber GmbH for inspection via your medical supply dealer.


6.1 Fault Displays on the Screen

Screen Display	Fault Description	Remedy
Error 20/ Error 21/ Error 22	The wheelchair is switched on but does not drive. Joystick defective.	Proceed as follows: <ul style="list-style-type: none"> • Switch off the Erivo (see Chapter 5.2). • Remove the battery pack (see Chapter 4.1). • Refit the battery pack (see Chapter 4.2). • Switch on the Erivo (see Chapter 5.1). • If the error code is still displayed: Contact your specialist dealer.
Error 24 / 42	The wheelchair is switched on but does not drive. Brakes are released.	Proceed as follows: <ul style="list-style-type: none"> • Switch off the Erivo (see Chapter 5.2). • Activate the parking brakes (see Chapter 5.11). • If the error code is still displayed: Contact your specialist dealer.
Error 40	The wheelchair is switched on but does not drive. Left parking brake deactivated.	Proceed as follows: <ul style="list-style-type: none"> • Activate the parking brake(s) (see Chapter 5.11). • If the error code is still displayed: Contact your specialist dealer.
Error 41	The wheelchair is switched on but does not drive. Right parking brake deactivated.	

Anzeige auf dem Bildschirm	Fehlerbeschreibung	Abhilfe
Error 80	The wheelchair is switched on but does not drive. Battery pack voltage too low.	<p>Proceed as follows:</p> <ul style="list-style-type: none"> • Switch off the Erivo (see Chapter 5.2). • Remove the battery pack (see Chapter 4.1). • Refit the battery pack (see Chapter 4.2). • Switch on the Erivo (see Chapter 5.1). • If the error code is still displayed: Contact your specialist dealer.
Error 81	The wheelchair is switched on but does not drive. Right motor defective.	
Error 82	The wheelchair is switched on but does not drive. Left motor defective.	
Error 83	The wheelchair is switched on but does not drive. Both motors defective.	
Error 84/ Error 85/ Error 86/ Error 87	The wheelchair is switched on but does not drive. An electronic switch in the controller between motors is defective.	
Error 88	The wheelchair is switched on but does not drive. Pre-charge error.	
Error 90	The wheelchair is switched on but does not drive. Battery pack voltage error.	
Error 91	The wheelchair is switched on but does not drive. Phase failure and disconnection error at left motor.	
Error 92	The wheelchair is switched on but does not drive. Phase failure and disconnection error at right motor.	
Error 93	The wheelchair is switched on but does not drive. Phase failure and disconnection error at both motors.	
Error 94/ Error 95/ Error 96	The wheelchair is switched on but does not drive. One or both motors blocked, causing overcurrent.	
E1	The wheelchair is switched on but does not drive. A control unit button was pressed for longer than 60 seconds.	
E6	The wheelchair is switched on but does not drive. Faulty communication between the control unit and the lower-level controller.	
E7	The wheelchair is switched on but does not drive. Invalid data in the configuration memory.	

6.2 Other Faults

Fault Description	Remedy
The Erivo cannot be switched on.	<ul style="list-style-type: none"> • Check whether battery pack [7] is correctly fitted and, if necessary, fit it correctly (see Chapter 4.2). • Check whether battery pack [7] is charged; charge if necessary (see Chapters 4.3 and 4.5 to 4.7). • If the battery pack is charged and correctly fitted: Contact your specialist dealer.
The Erivo does not drive straight.	<ul style="list-style-type: none"> • Contact your specialist dealer.
The Erivo suddenly slows down at high speed.	
One or both Erivo motors make unusual noises.	

7. Automatic Self-Shutdown

7.1 Self-Shutdown Due to Internal Faults

In the event of a system malfunction or overloading of the circuit (e.g. caused by high-load driving on inclines), the Erivo automatically switches off. To recommission it, restart the system (see Chapter 5.1):

Press the On/Off / Anti-theft lock button [20] on the control unit [2].

7.2 Self-Shutdown During Non-Use / Cancelling Self-Shutdown

If the Erivo is switched on but not used, it automatically switches itself off after 30 minutes.


- Press the On/Off / Anti-theft lock button [20] to cancel the self-shutdown and restart the system


8. Checking Tyre Pressure

The electrically driven rear wheels of the Erivo are fitted with pneumatic tyres.

- Regularly check the tyre pressure of both tyres (see the table in Chapter 17) and top up via valve [65] if necessary.

The two castors (front wheels) are puncture-proof and therefore do not require air.

9. Transferring

- Switch off the Erivo.
This automatically engages the parking brakes, preventing unintentional movement.
- If necessary, swing one or both armrests [3] upward for transferring.
- Lean on one or both armrests [3], but never on the control unit.


Before transferring, the Erivo must be switched off and thus the parking brakes engaged, to prevent unintentional movement of the wheelchair.

10. Removing / Placing the Seat Cushion

The seat cushion [14] is connected to the wheelchair seat by two hook-and-loop straps [68].

- Remove the seat cushion [14] if you wish to clean it, or place it back on the seat after cleaning.

11. Transport

Handle the Erivo carefully during transport, do not drag it along the floor, do not drop it, and avoid impacts etc.

11.1 Folding and Unfolding the Wheelchair


Safety notice — Crush hazard!

When folding or unfolding, hold the wheelchair only at the positions specified in these chapters. Never reach into the area of mechanical parts that move together when folding or unfolding, and never touch the pinch points marked on the illustrations and/or on the wheelchair.

Areas where there is a crush hazard are marked on the wheelchair with stickers bearing the symbols shown here.

11.1.1 Folding the Wheelchair

For transport purposes, an accompanying person can fold your Erivo.

- Switch off the Erivo, see Chapter 5.2.
- Remove the battery pack for transport to prevent the Erivo from accidentally switching on and moving, see Chapter 4.1.
- Place the Erivo on a level surface - never on sloping terrain.
- Empty the storage basket [11] before folding to avoid damaging items inside.
- Stow the safety belt as follows:

- Open the safety belt.
- Roll both parts together.
- Secure each part with the rubber band on each part (see illustration).


Alternatively, you can stow the safety belt as follows:

- Close the safety belt.
- Lay the safety belt on the seat cushion [14].
- Stand in front of the Erivo facing the seat cushion [14] and backrest.
- Press the lever [47] of the folding mechanism fully against the Erivo frame and **hold it firmly**.
- Simultaneously pull the push handle [5] at the upper end of the backrest towards you and downward. Hold one armrest [3] to prevent the Erivo from rolling during folding.


The lever [47] can be released after you have pulled the push handle [5] slightly towards you.

The Erivo is fully folded when you hear a quiet „click“ and the wheelchair can no longer be pulled apart without operating rod [48].

- Fold the footrest [13] upward so that it rests close against the push handle [5]. Only fold the footrest up after the Erivo is already folded (see the sticker on the underside of the footrest [13]).


NOTE: WHEN FOLDING THE WHEELCHAIR PLEASE FOLD THE FOOTRESTS LAST.

11.1.2 Unfolding the Wheelchair

To return the folded Erivo to operational status, an accompanying person can unfold it:

- Place the folded Erivo on a level surface - never on sloping terrain.
- Fold the footrest [13] downward.
- Stand in front of the Erivo facing the seat cushion [14] and backrest.
- Press the rod [48] firmly downward until you hear a quiet „click“.
- Pull the push handle [5] slightly upward.
The rod [48] can be released once the push handle [5] has been pulled slightly upward.
- Pull the push handle [5] further upward and simultaneously press the seat surface [14] downward.
The Erivo is fully unfolded when you hear a quiet „click“ and the wheelchair stands on all four wheels.
- Before commissioning, verify that the Erivo is fully unfolded and can no longer be folded without operating lever [47] (see Chapters 11.1 and 11.1.1).

11.2 Safety Instructions for Loading the Wheelchair

Observe the safety and hazard notices in Chapters 3.4 to 3.6

11.2.1 Loading Without a Ramp or Lifting Device - Folded

- Preparation: Read and observe the notes in Chapter 11.1.1.
- Ensure the Erivo is completely folded and can no longer be pulled apart without operating rod [48] (see Chapter 11.1.2).
- Load the Erivo while folded with the help of an accompanying person.
- When loading or unloading the Erivo, never hold it by both armrests [3] as these cannot be locked (crush hazard!).
- Never hold the Erivo by parts of the footrest [13] as this cannot be locked (crush hazard!).

11.2.2 Loading via Ramp or Lifting Device - Unfolded - by Accompanying Person

- Switch off the Erivo (see Chapter 5.2).
- Activate Push Mode of the wheels (see Chapter 5.10).
- Maintain sufficient clearance from the edges of the ramp or lifting device when pushing.
- Use only ramps or lifting devices with appropriate load capacity.
- Only use a ramp or lifting device if you are fully familiar with its operation according to the manufacturer’s specifications.

11.2.3 Loading via Ramp or Lifting Device - Unfolded - by Wheelchair User

- Navigate a ramp or lifting device with extreme caution and exclusively at minimum speed.
- On a lifting device: maintain sufficient clearance from the edge.
- Do not drive curves on a ramp.
- Use only ramps or lifting devices with appropriate load capacity.
- Only use a ramp or lifting device if you are fully familiar with its operation according to the manufacturer's/distributor's specifications.

11.3 Transporting the Wheelchair in a Vehicle

The Erivo and its components (battery pack, control unit) must be secured during transport in a vehicle so that parts cannot become a hazard to vehicle occupants in the event of braking. Components must never be stored unsecured! If national regulations of your country apply regarding securing the Erivo or its components, these take precedence and must be observed. If national regulations of your country apply regarding securing the Erivo or its components, these take precedence and must be observed.

The Erivo must not under any circumstances be used as a vehicle seat.

12. Storage, Care, Recommissioning, Maintenance and Disposal

12.1 Storing the Wheelchair

Observe the following when storing your Erivo for an extended period (e.g. several months):

- Protect the Erivo components with plastic sheeting to keep out moisture.
- Store the Erivo and its components in a dry room.
- Ensure no moisture enters or condenses on the Erivo and its components.
- Ensure the Erivo and its components are not exposed to prolonged sunlight (e.g. through windows).
- Ensure unauthorised persons, in particular children, do not have unsupervised access to the storage room.
- Observe the notes on storing the battery pack [7] in Chapter 4.8.
- If folding is required for storage: Ensure nothing is trapped during folding and that the storage basket [11] is empty.
- If storing the Erivo folded: Store the electric wheelchair upright leaning against a wall to prevent deformation of the frame and tyre flat spots (see illustration).
- Clean the Erivo and its components before recommissioning. In particular, no metallic particles should be present in or on the plugs [52] and [53] of the charger, the socket on the battery pack [49], or the socket [49] of the control unit.
- Check whether maintenance (see Chapter 12.4) is required before recommissioning and arrange for it accordingly.

12.2 Caring for the Wheelchair

Ingress of water can destroy the drive.

Therefore, never clean the Erivo or its individual components with running water, a high-pressure cleaner, or in a car wash.

Always ensure that neither liquids nor moisture enter the components, especially the wheels and battery pack!

Clean the battery pack exclusively with damp (not wet!) cloths. Under no circumstances must liquid enter the battery pack (fire risk!).

If chemical splashes land on the Erivo: Immediately rinse the area carefully with clean water. Dry with a soft cloth. If necessary, contact an Alber representative or an Alber-authorized specialist dealer.

When cleaning, ensure no water enters Erivo components such as the battery pack or wheels.

Alber GmbH disclaims all liability for damages or consequential damages caused by water ingress. Such cases are also not covered by warranty.

Regular cleaning

- Switch off the Erivo before cleaning or maintenance work to prevent unintentional movement.
- Regularly check the Erivo for cleanliness, especially after driving over grass, sand, rough surfaces and potholes, and after intensive use.
- Do not expose the Erivo to heat, cold or moisture for extended periods.
- Keep the control unit always clean.
- Use only slightly damp cloths for all cleaning; never abrasive agents or aggressive cleaning products such as steel wool.
- Never apply liquid cleaners or solvents to the control unit, battery pack, charger or electrical connections.
- When cleaning, pay particular attention to the joystick. It must always return to its central position. The rubber surround must also be undamaged to prevent moisture ingress.
- Regularly clean the wheels and axles.
- If you notice that parts of the Erivo are no longer fully functional, screws are missing, exposed wires are visible, or other externally visible damage (e.g. scratches, deeper frame damage, delamination, deformation) is apparent, contact your specialist dealer immediately for a function check or repair.

Annual cleaning

The following components should be thoroughly cleaned at least once annually (or as required).

- Armrest: Wipe with a damp cloth and gentle cleaning agent.
- Seat cushion and backrest: First use a vacuum cleaner, then a damp cloth with gentle cleaning agent. Never use abrasive cleaners as these can damage the seat cushion.
- Control unit: Wipe with a damp cloth, then dry with a dry cloth. Never use abrasive cleaners.

12.3 Recommissioning

If the Erivo was provided to you by your health insurance and you no longer need it, please contact your health insurance, an Alber representative, or your medical supply dealer. Your Erivo can then be easily and economically recommissioned.

Before each recommissioning, the Erivo should be serviced.

In addition to the cleaning instructions in Chapter 24.2, before recommissioning, all externally accessible plastic parts of the Erivo components must be disinfected. Use only alcohol-based wipe disinfectants approved by the Robert Koch Institute (RKI); follow the manufacturer's recommended contact time and concentration.

Example: Bacillol AF, contact time 15 minutes.

12.4 Maintenance

The legislator has regulated the maintenance of medical devices in the Medical Device Operator Regulation (MPBetreibV) § 7. Maintenance measures, especially inspections and servicing, are required to continuously ensure the safe and proper operation of medical devices.

For the maintenance of our products, an interval of two years has proven useful based on market observation under normal operating conditions. This guideline of two years may vary depending on the usage level and user behaviour. Assessment of usage level and user behaviour is the responsibility of the operator. In any case, we recommend clarifying cost coverage for maintenance work with the relevant payers (health insurers) in advance, especially with regard to any existing contracts.

12.5 Disposal

This device, its battery pack and accessories are long-lasting products; however, they will eventually reach the end of their service life and must be disposed of properly.

They may contain substances classified as harmful to the environment under current EU regulations if disposed of at locations not designated for this purpose under current national regulations.

The „crossed-out wheeled bin“ symbol (in accordance with the WEEE Directive and EU Battery Regulation) is on this product as a reminder of your obligation for proper and compliant disposal.

They also contain (sometimes rare) valuable materials that can be recycled and returned to the material cycle.

Please act in an environmentally responsible manner and take this product to your local recycling facility at the end of its service life — usually municipal recycling centres.

Please inform yourself about the regulations applicable in your country regarding disposal.

Alternatively, the Erivo can be returned to Alber or Alber specialist dealers for proper and environmentally friendly disposal.

13. Product Safety Notices

All components of your Erivo have been subjected to numerous functional tests and extensive checks. Should unforeseeable impairments in the operation of the Erivo nonetheless occur, relevant safety information for customers would be available promptly on the following websites:

- Website of Alber GmbH (distributor of the Erivo)
- Website of the Federal Institute for Drugs and Medical Devices (BfArM)

In the event of any necessary corrective measures, Alber will inform the medical supply trade, which will then contact you.

14. Product Lifespan

Regarding the estimated lifespan of this product, we assume an average of five years, provided the product is used within its intended use and all maintenance and service requirements are met. This lifespan can be exceeded if the product is carefully handled, maintained, cared for and used, and no technical limitations arise from advances in science and technology. However, the lifespan may be significantly shortened by extreme use and improper handling. The normatively required specification of lifespan does not constitute an additional guarantee.

15. Warranty, Guarantee and Liability

15.1 Defect Warranty

Alber warrants that the Erivo is free of defects at the time of delivery. These warranty claims become statute-barred 24 months after delivery of the Erivo.

15.2 Durability Guarantee

Alber provides a 24-month durability guarantee on the Erivo.

The durability guarantee does not cover

- Devices whose serial numbers have been altered, defaced or removed.
- Wear parts such as tyres and controls
- Defects due to natural wear, mishandling, in particular defects due to non-observance of these instructions for use, accidents, negligent damage, fire, water damage, force majeure and other causes outside Alber's area of influence.
- Maintenance work required by daily use (e.g. tyre replacement).
- Device inspection without defect finding.

15.3 Liability

Alber GmbH, as the distributor of the Erivo, is not responsible for its safety if:

- the Erivo is improperly handled
- the Erivo is not serviced every 2 years by an authorised specialist dealer or Alber GmbH
- the Erivo is commissioned contrary to the instructions of these instructions for use
- the Erivo is operated with insufficient battery charge
- repairs or other work are carried out on the Erivo by unauthorised persons
- foreign parts are fitted or connected to the Erivo
- parts of the Erivo are removed

16. Important Legal Notice for Users of this Product

Incidents occurring as a result of product malfunctions that lead to serious personal injury must be reported to the manufacturer and the competent authority of the EU Member State in which the user is established.

In the Federal Republic of Germany these are

- MedNet EC-REP GmbH
Borkstrasse 10
48163 Muenster
- Federal Institute for Drugs and Medical Devices (BfArM),
Kurt-Georg-Kiesinger-Allee 3
53175 Bonn
www.bfarm.de

In EU Member States, notify

- Your Alber representative (address on reverse of this IFU)
- The authority responsible for incidents involving medical devices in your country

An overview of the competent authorities can be found online at: <http://ec.europa.eu/growth/sectors/medical-devices/contacts/>

17. Technical Data

Wheelchair

Range (1):	Up to 13.4 km with standard battery pack 12 Ah (Application Class A per EN 12184) Up to 21.5 km with optional battery pack 20 Ah (Application Class A per EN 12184)
Rated gradient:	Up to 9° (16%)
Obstacles:	Maximum permitted obstacle height/depth: 4 cm
Speed:	1,3 - 6 km/h Applicable national road traffic laws must be observed!
Motor rated power:	2 x 250 W
Operating voltage:	24 V
Operating temperature:	-10°C to +50°C
Permitted humidity:	Relative humidity: 20 % to 95 %
Permitted atmospheric pressure:	70 kPa to 110 kPa
User weight:	Maximum 150 kg

Battery Pack

	Standard Battery Pack	Optional Battery Pack
Cell type:	Lithium-ion cells, leak-proof and maintenance-free	
Operating nominal voltage:	29.4 V	29.4 V
Nominal capacity:	12 Ah	20 Ah
Nominal energy:	288 Wh	480 Wh
Charging temperature:	0°C to 45°C	0°C to 45°C
Operating temperature:	-20°C to +60°C	-20°C to +60°C
Ingress protection:	IPX5	IPX5

System Storage

Storage temperature:	-20°C to +60°C
Relative humidity during storage:	max. 95%
Atmospheric pressure during storage:	56 to 110 kPa

Weights of Individual Parts

Wheelchair without battery pack and without bag:	16.4 kg	
Bag:	1.0 kg	
Battery pack:	2.5 kg (standard)	3.4 kg (optional)
Total weight:	19.9 kg	20.8 kg

Wheel Tyres

Designation	Type	Size (inches)	Tyre pressure in bar / kPa / psi
Drive wheels: Xuan Cheng	10 x 2.125 HY-202-10	10"	3.1 bar / 310 kPa / 45 psi
Castors: PU front wheel tyres	not specified	8"	Puncture-proof


(1) Range varies depending on terrain and prevailing driving conditions. Under optimal conditions (flat terrain, fully charged battery, ambient temperature of 20 °C, steady driving etc.) the stated range can be achieved. At lower temperatures the range may be less than indicated above.

Ladegerät

Model:	YS90-2943000
Mains voltage:	100 - 240 V AC 50/60 Hz
Output power:	88.2 W
Output voltage:	29.4 V DC
Output current:	3.0 A
Ingress protection:	IPX0
Protection class:	n.a.
Ambient temperature:	Operation: -20° C bis +50° C Storage: -40° C to +80° C
Humidity:	Operation: n.a. Storage: n.a.
Air pressure:	Operation: n.a. Storage: n.a.
Charger plug polarity:	

Block diagram

CE The Erivo fulfils all applicable chapters of EN 12184 for electrically powered wheelchairs, as well as the relevant normative references of this standard, in particular ISO 7176-14 and EN 1021-2 (or equivalent ISO 8191-2), and thus complies with EU Medical Device Regulation 2017/745 (MDR) and is classified as a Class I medical device under this regulation.

Subject to technical and design changes due to continuous development.

18. Labels and Markings

Labels are located on the Erivo components and contain various product information or warnings. The following illustrations are examples only!

Electrically Powered Wheelchair ROBOTEER E80	Type classes: Class A Max. mass of the user: 150Kg Input: 24V 3A	CE # IPX5 LOT MD SN	BBR - E80-01 J2025112030 B085BX194 2025-12-13
Jiangsu Bangbang Intelligent Technology Co., Ltd. Building 1, No.2 Hongren Road, Tianning District, Changzhou Jiangsu, China Email: info@roboteer.com	EU REP SUNGO Europe B.V. Fascinatio Boulevard 522, Unit 1.7.2909 VA Capelle aan den IJssel Netherlands ec.rep@sungoglobal.com/ yan.zhang@sungoglobal.com	UDI	

Manufacturer System Label

The system label is located on the right-hand side of the frame, beneath the seat. The system label provides key technical specifications for the Erivo.

Systemnr./System No. R202699999
 Produkt/Modell [Product/Model] erivo carbon R20
 Artikelnr. [Item no.] 1596014

Distributor System Label

The system label is located on the vertical (rear) section of the battery compartment. The system label enables the system to be traced if necessary.

WARNING **WARNING**
 Hands pinching Hands pinching

Crush Hazard Warning Label.

The labels are located on the frame, above the levers [47] used to fold the wheelchair. There is a risk of crushing your hands at the marked points on the Erivo, particularly when folding the wheelchair.

PERMANENT MAGNET DC BRUSHLESS MOTOR

MODLE	BBR-M-04	POWER	250W
VOLTS	24VDC	REDUCTION RATIO	17.2
R. P. M	167	INS. F	IP55

Jiangsu Bangbang
 Intelligent Technology Co., Ltd.

P/N: 016. 004. 0013664
 P/O: DXC-LCIDZBU250701196
 S/N: 20250809L0142

Motor Data Label

The label is located on both motors. It provides information about their technical specifications.

NOTE: WHEN FOLDING THE WHEELCHAIR PLEASE FOLD THE FOOTRESTS LAST.

Folding Sequence Label - Footrest


The label is located on the underside of the footrest. It advises you that when folding the Erivo, you must follow the correct sequence for folding and raising the footrest; see Chapter 11.1.1.

WARNING

 It can not be fixed on the vehicle as a seat

Warning Label

The label is located on the frame above the left rear wheel. It informs you that you must not use the Erivo as a vehicle seat.


Parking Brake Activation / Release Label

The labels are located on the cover of the electric motors, above the parking brakes [8]. They indicate that the parking brakes [8] are engaged when in the upper position and disengaged when in the lower position.

Battery Label

The label is located on the battery pack. It contains technical and other information, as well as warnings regarding the handling of the battery pack.

Battery Lock / Unlock Marking

One marking (1) is located on the left-hand side of the battery pack [7], the other marking (2) on the right-hand side of the battery pack [7]. The markings indicate the position in which the buttons [45] are locked or unlocked.

Charger Label

The label on the charger provides key technical specifications and instructions on how to use the charger.

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Manufacturer:
SHANGHAI BANGBANG ROBOTICS CO. LTD.
District Shanghai
Zhongchen Road, Songjiang
Room 501, Building 3, No. 188
China

UK	RP
----	----

Responsible Person United Kingdom:
Medimap Ltd
2 The Drift
Suffolk
Thurston IP31 3RT
United Kingdom

Appendix

B - Inspection Checks

Check the functionality of the Erivo before every journey. If malfunctions occur, please contact your authorised specialist dealer.

Inspection Point	Check Step	Remedy if Malfunction
Control unit	Check that the battery charge displays [30], [31] illuminate on the screen after switch-on.	Is the battery charged and inserted?
	Check for error messages on screen [23] (see Chapter 18).	Rectify faults as per this IFU.
	Check battery pack capacity (see Chapter 4.4).	Fully charge battery pack before starting.
	Check correct function of On/Off / Anti-theft lock button [20] and joystick [24].	Rectify any faults as per this IFU.
	Check joystick position (with control unit switched off!). Ensure it is not bent/damaged and returns to centre when released. Visually inspect rubber boot at joystick base. Do not attempt to repair the joystick.	Contact your specialist dealer if damage/malfunction found.
	Check function of horn button [22].	Contact your specialist dealer.
Motors	Check motors for unusual noises, overheating, malfunctions.	Contact your specialist dealer.
Brakes	In electric motor mode, check that brakes hold the Erivo when joystick is released.	Contact your specialist dealer.

Inspection Point	Check Step	Remedy if Malfunction
Operating mode	Check that motors switch smoothly between electric motor mode and Push Mode.	Contact your specialist dealer.
Mechanics	Check that folding/unfolding works correctly.	Contact your specialist dealer.
	Check whether frame, armrest and footrest connections have deformed or loosened.	Contact your specialist dealer.
	Check all parts are properly bolted and no parts (e.g. screws) are missing.	Contact your specialist dealer.
	Check removable parts are correctly connected and secured.	Rectify faults as per this IFU.
Battery pack	Charging at control unit: Check that display on control unit and LED indicator on charger illuminate after connecting battery.	Check cable connection from charger to control unit. If unresolved, contact your specialist dealer.
	Charging at battery pack: Check that LED indicator on charger illuminates after connecting battery.	Check cable connection from charger to battery. If unresolved, contact your specialist dealer.
	The charger LED does not illuminate green when the battery is fully charged.	Contact your specialist dealer.
Wheels & tyres	Check tread depth on all wheels.	Tyres with tread depth below 1 mm must be replaced. Contact your specialist dealer. The Erivo must not be operated with worn tyres (no visible tread).
	Check tyres for cracks or deformations. Remove any foreign objects from tyres.	If cracks/deformations found, tyre must be replaced. Contact your specialist dealer.
	Check tyre pressure.	Is air constantly escaping? Contact your specialist dealer.
	Check that wheels rotate without wobbling.	Contact your specialist dealer.
Charger	Check charger for correct function, especially LED indicators.	If malfunctioning, charger must no longer be used. Contact your specialist dealer.
Cables	Visually inspect any exposed cables. Ensure they are not frayed, cut, or have exposed wires.	Contact your specialist dealer if any wiring harness issue is found.

If faults cannot be resolved, contact your specialist dealer immediately.

Ihre Alber-Vertretung / Your Alber representative / Votre représentation Alber /
Vostra rappresentanza Alber / Su representación Alber / Din Alber representant /
Din Alber-agenturene / Uw distributeur Alber / Deres Alber-repræsentation

**SHANGHAI BANGBANG ROBO-
TICS CO. LTD.**

District Shanghai
Zhongchen Road, Songjiang
Room 501, Building 3, No. 188
China

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MedNet SWISS GmbH